

MACROMAN™

ΔΙΑ•М
современная лаборатория

www.dia-m.ru
заказ on-line



*Unique Comfort
& Controlled Pipetting*



Stable, safe & easy to pick up

- > Integrated **bench stand** adjustable for **left and right** handed users.
- > **Keep your MACROMAN handy** and ready for use, without risk of contamination.

Soft & fatigue free pipetting

- > **Soft-touch material** for maximum finger comfort.
- > **Large dosing button** to decrease pushing efforts.



MACROMAN™



Optimal comfort

Integrated, adjustable stand for easy handling when pipette attached

Large & soft-touch material dosing button

Light weight

Quiet operation & **vibration free**

Ergonomic-shape handle with grip

Safety & easy maintenance

All liquid contact parts are autoclavable

Easy-to-replace filter

Replaceable batteries

Performance

For all plastic & glass pipettes from 1 to 100mL

Continuously adjustable speed: 25mL in 4sec

Environmentally friendly NiMH battery
for **8 hours operation**

Precise & smooth dosing button control: drop-by-drop, gravity mode
and motorized speed.

Convenient accessories

Horizontal or vertical & magnetic charging stand (with complete mounting kit) for busy laminar flow cabinet.

Unique adaptor for D5000 and D10mL DIAMOND® tips.



Charging stand



MACROMAN

Weight	194g
Battery	2x 1.2V AAA NiMH, removable, 8-hour operation, rechargeable during use, discharged & charging light.
Autoclavability	Nose piece, Silicone pipette mount & filter rubber.
Filter	Replaceable, 0.45µm or 0.20µm standard & sterile.
Pipette types	Glass or plastic from 1mL to 100mL
Power supply	100-240V / 50-60Hz/ EU, US, UK, JAP or AUS plug version
Certification	CE, EMC, RoHS, WEEE
Contents	Complete MACROMAN + holder + 1 extra 0.45µm sterile filter + identification labels + EU, US, UK, JAP or AUS plug version

Description	Reference Number
MACROMAN, EU plug	F110751
MACROMAN, US plug	F110752
MACROMAN, UK plug	F110753
MACROMAN, JAP plug	F110754
MACROMAN, AUS plug	F110755
Charging stand	F1077501
0.45µm standard filter	F1077504
0.45µm sterile filter	F1077525
0.20µm standard filter	F1077503
0.20µm sterile filter	F1077526
Nose piece	F1077506
Filter rubber	F1077502
Pipette mount	F1077505
Adapter for D5000 and D10mL DIAMOND® Tips	F161287

Gilson provides 1-year warranty after the date of purchase

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

